



## Drug of Abuse Oral Fluid Controls

- Synthetic oral fluid matrix simulates patient samples
- Gravimetrically manufactured to target levels
- Stable for up to 2 years frozen (unopened) with 30-day open stability
- Ready to use liquid controls in convenient 5 or 20 mL size vials
- Designed to evaluate the performance by all oral fluid testing methodologies



We also offer a complete line of controls for Urine Drugs of Abuse (DOA) testing

Vendor Catalog#	Liquid Oral Fluid Quality Controls	Kit Size
702030	<b>SALIVABUSE ORAL SKREEN, POSITIVE &amp; NEGATIVE</b> DOA 15 PANEL Drug Control	2 X 20 mL
20200101-6PK	<b>SALIVABUSE ORAL SKREEN, NEGATIVE</b>	6 X 20 mL
20200109-6PK	<b>SALIVABUSE ORAL SKREEN, POSITIVE</b> DOA 15 PANEL Drug Control	6 X 20 mL
20200100	<b>SALIVABUSE ORAL SKREEN, NEGATIVE</b>	500 mL
20200200	<b>SALIVABUSE ORAL SKREEN, NEGATIVE with STABILIZER</b>	500 mL
18100156	<b>SALIVABUSE LIQUID CONTROL, AUNZS19*</b> CUTOFF -50%	20 mL
18100157	<b>SALIVABUSE LIQUID CONTROL, AUNZS19*</b> CUTOFF +50%	20 mL

\*Formulation set to Australia and New Zealand 2019 regulations

Contact your Biochemical Diagnostics Sales Rep at [KOVA.SUPPORT@LGCGROUP.COM](mailto:KOVA.SUPPORT@LGCGROUP.COM)  
for additional product information, visit our website at [kovaintl.com](http://kovaintl.com)



## SALIVABUSE® Liquid Control Oral Fluid, ORAL SKREEN

Target Values (ng/mL)			
Test	Constituents	NEGATIVE	3X POSITIVE
THC	Delta-9-THC (Parent THC)	0	120
COC	Cocaine (Parent)	0	90
COC	Benzoyllecgonine	0	90
PCP	Phencyclidine (PCP)	0	30
OPI	Morphine	0	120
OXY	Oxycodone	0	90
AMP	Amphetamine	0	150
MDMA	MDMA	0	150
MET	Methamphetamine	0	150
BAR	Secobarbital	0	150
BZO	Oxazepam	0	150
MTD	Methadone	0	90
COT	Cotinine	0	90
ALC	Ethanol (mg/dL)	0	75
BUP	Buprenorphine	0	30

## SALIVABUSE® Liquid Control Oral Fluid, AUNZS19

Target Values (ng/mL)			
	18100156 (20mL)		18100157 (20mL)
CONSTITUENTS	CUTOFF -50%	CUTOFF	CUTOFF +50%
Oxazepam	5	10	15
d-AMPHETAMINE	25	50	75
d-METHAMPHETAMINE	25	50	75
$\Delta^9$ -tetrahydrocannabinol (Delta-9-THC, parent drug)	7.5	15	22.5
Benzoyllecgonine	25	50	75
Morphine	25	50	75
Oxycodone	20	40	60

